


Manza Heads

TECHNICAL FEATURES



Manza Heads

TECHNICAL FEATURES



All heads are supplied in with a 12/14 taper specification per Ceramtec.
 Offsets are identical to Ceramtec to allow identical interchangeability between bearing surfaces.
 All Manza heads are a skirtless design to maximise Range Of Motion (ROM).
 Heads available in either CoCr, High Nitrogen Stainless Steel or BioloX Delta.

Head Diameters (mm)	Offset				Material		
	S	M	L	XL	CrCo	Stainless Steel	BioloX Delta
28	S	M	L		CrCo	Stainless Steel	BioloX Delta
32	S	M	L	XL	CrCo	Stainless Steel	BioloX Delta
36	S	M	L	XL	CrCo	Stainless Steel	BioloX Delta
40	S	M	L	XL	CrCo	Stainless Steel	BioloX Delta

Proven, Unique, Enduring

Handwriting practice area with 20 horizontal dashed lines.

